



neoLab® Thermo-/Hygrometer, Max./Min. function -10°C/+60°C

Product Images



Description

Electronic thermometer and hygrometer with maximum/minimum function, for standing or hanging.

Technical data

Temperature: -10°C to +60°C

Rel. Humidity: 10-99%

Resolution: 0.1°C

Accuracy: $\pm 1^\circ\text{C}$, $\pm 5\%$

Display: temperature and humidity simultaneously, switchable °C/°F

Memory function: maximum/minimum recall for temperature and humidity

Power supply: 2 micro batteries 1.5 V incl.

Dimensions: 110 x 70 x 20 mm

Weight: 113 g

Additional Information

No.	2-5437
Manufacturer (Brand)	neoLab
EAN	4058072004682
Color	Weiß
Temp MAX	60 °C
Temp MIN	-10 °C
Wide	70 mm
Height	110 mm
Length	70 mm
Depth	70 mm
Display (reading, setting)	digital
Weight	74 g
Counting and measuring devices Features	mit interner Kalibrierung inklusive Batterien/Akku
(Combined) functions	Thermometer Hygrometer
Resolution °C	0.1 °C
Energy supply	Batterie
Resolution rH %	1.00 %
Calibration type	automatisch
measured size(s)	Temperatur Luftfeuchte
Measuring unit	°C °F rF %
Measuring accuracy rH %	5.00 rF
Measuring accuracy °C	1 °C
rF MAX	99.00 rF
rF MIN	10.00 rF

