



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Hand piece counter 0-9999, with ring

€74.30
plus VAT &
Shipping

Product Images



Description

Hand counter with retaining ring, for counting pieces up to 9,999.

Mechanical, very reliable counter.

White digits, 4-digit, on black background.

Rotating zero position with knurled head.

Housing made of zinc die-cast alloy.

Operating temperature: 0 to +60°C

Housing diameter : 47 mm

Weight: approx. 115 g

Quality made in Germany

Additional Information

No.	1-7509
Manufacturer (Brand)	neoLab
EAN	4058072046460
diameter	47 mm
Color	silver
Material	Steelgalvanized
Display (reading, setting)	digital
Operation via	Knobs/Buttons
Weight	115 g
Number of digits	9999
Functionality	mechanical

