



neoLab® Mechanical counter, 8-fold

Product Images



Description

Mechanical counter, especially suitable for counting cells, particles, etc.

Stickers for blood cells are included.

Each key counts up to 999, separate display for each key.

Another viewing window shows the total number of all counts.

A bell sounds each time a hundredth step is reached.

Order No. 2-5097
8-fold counter
Dimensions: 320 x 76 x 56 mm (L x W x H)

Additional Information

No.	2-5097
Manufacturer (Brand)	neoLab
EAN	4058072004347
Color	anthrazit
Wide	76 mm
Height	56 mm
Length	320 mm
Depth	320 mm
Display (reading, setting)	analog
Operation via	Drehknöpfe/Tasten
Counting and measuring devices Features	mit Fuß
Number of digits	999
Number of counters	8
Functionality	mechanisch

