



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific LF Pocket Tester ECTestr 11, 0-2000 μ s, 0-20 mS, AutoCal

€158.00
plus VAT &
Shipping

Product Images

Description

This series covers the entire measuring range from pure water to waste water. The testers have a 1% measuring accuracy over the entire measuring range. They can be set to the measurements with up to 3 calibration points, which allows maximum accuracy to be achieved.

Order No. 8-1305

ECTestr 11 with AutoCal

Measuring range: 0-2000 μ S/ 0-20mS

Additional Information

No.	LD-3965
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	8-1305
EAN	4066292033726
Color	Yellow White
Temp MAX	50 °C
Temp MIN	0 °C
Wide	65 mm
Height	50 mm
Length	185 mm
Depth	185 mm
Display (reading, setting)	digital
Weight	200 g
for ambient temperature	0 to 50 °C
Counting and measuring devices Features	waterproof with automatic temperature compensation (ATC) batteries/battery included
Resolution °C	0.1 °C
(Combined) functions	Conductivity meter
Energy supply	Battery
Conductivity MAX	20000.000 µS/cm
Conductivity MIN	0.000 µS/cm
Resolution µS/cm	0.100 µS/cm
Calibration type	automatic
measured size(s)	Conductivity
Spare part	LD-3964
Measuring unit	µS/cm
Measuring accuracy µS/cm	1.00 %

