



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Mettler-Toledo pH meter FiveEasy , F20 standard kit

€740.00
plus VAT &
Shipping

Product Images



Description

FiveEasy™ pH F20 standard kit

Compact entry-level instrument for standard pH and mV/redox measurements with unbeatable price/performance ratio.

Applicable in many fields (e.g. water/environmental analysis, food and beverage industry, education and others).

-3-point calibration, 5 predefined buffer groups

- Large LCD display (everything at a glance), universal power supply

Including pH electrode LE438

Range of application:

Suitable for a wide range of aqueous samples

Additional Information

No.	TD-0164
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105006626
EAN	4058072983017
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	whitelight gray
Material	Acrylonitrile butadiene styrene (ABS)
Temp MAX	100 °C
Temp MIN	0 °C
Wide	227 mm
Height	70 mm
Depth	147 mm
Display (reading, setting)	Display 4.3" LCD Color
Weight	0.63 kg
Type interface	USB RS232
for ambient temperature	0 to 40 °C
Counting and measuring devices Features	with internal calibration with integrated data memory with automatic temperature compensation (ATC) with automatic buffer recognition with sensor with interface with limit value setting with alarm function including batteries/battery pack
Resolution °C	0.1 °C
(Combined) functions	pH meter thermometer
Energy supply	Power grid
Calibration type	at the touch of a button
pH MAX	14.0 pH
pH MIN	0.0 pH
Data transmission	Cable
Shape connection	straight

Manufacturer	Mettler-Toledo OnLine
Measuring accuracy °C	0.1 °C
Type sensor	PH electrode LE438
Type plug	BNC Cinch

