



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Mettler-Toledo oxygen meter Seven2Go polarographic DO Meter S4 standard kit

€1,100.00
plus VAT &
Shipping

Product Images



Description

Seven2Go portable oxygen meter (DO meter) S4-Meter

Intuitive routine instrument (IP67) for mobile measurements of dissolved oxygen (DO) in laboratory, production and field with a wide range of functions:

- very simple, intuitive operation, robust and ergonomic design
- large multifunction display: everything at a glance
- ISM® (intelligent sensor management)
- data logger for 200 data sets
- automatic air pressure compensation
- automatic, timer-controlled or manual measured value acquisition
- water and dust protected (IP67)
- operation: AA batteries or rechargeable batteries

Including InLab® Electrode 605-ISM

Additional Information

No.	TD-0132
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105002652
EAN	4066292341043
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	whitelight gray
Material	Acrylonitrile butadiene styrene (ABS)
Temp MAX	60 °C
Temp MIN	0 °C
Wide	70 mm
Height	35 mm
Depth	222 mm
Display (reading, setting)	Display 4.3" LCD Color
Weight	0.27 kg
for ambient temperature	0 to 40 °C
Counting and measuring devices Features	waterproof with internal calibration with integrated data memory with automatic temperature compensation (ATC) with sensor with limit setting with alarm function batteries/battery included
Resolution °C	0.1 °C
(Combined) functions	Manometer Thermometer Oxygen meter Saturation meter Oxygen
Pressure MAX	1100.000 mbar abs.
Pressure MIN	500.0000 mbar abs.
Energy supply	Battery Battery
Calibration type	at the touch of a button
measured size(s)	Temperature Pressure Dissolved oxygen Oxygen saturation
Shape connection	straight
Manufacturer	Mettler-Toledo OnLine

Measuring unit	mbar abs. mg/l ppm °C
Measuring accuracy mbar	10 mbar
Measuring accuracy °C	0.2 °C
Type sensor	InLab 605-ISM
Type plug	BNC Cinch

