



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle flow indicator made of SAN f. hoses 6-11 mm inner Ø

€11.15 As low as
€10.25
plus VAT &
Shipping

Product Images



Description

Flow indicator for accurate visual indication of flow velocities.

It works in both flow directions and in any installation position.

The impeller rotates from approx. 150 ml/min.

Dimensions: 88 x 40 x 15 mm (L x W x D)

Body Ø: 41 mm

Additional Information

No.	KL-1544
Manufacturer (Brand)	Bürkle
old neoLab article no.	1-1340
EAN	4058072274344
Transport temperature	Room temperature
diameter	41 mm
Color	transparent (clear)
Material	Styrene-acrylonitrile (SAN)
Wide	40 mm
Length	88 mm
TBST MAX	30 °C
MAX operating pressure	2 bar abs.
for DM hose inside 1	6.0 to 11.0 mm

