




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Amarell general purpose thermometer, enclosed type, -35+50:1°C, capillary prismatic

**€24.20** As low  
as  
**€15.95**  
plus VAT &  
Shipping



## Description

---

General purpose thermometer, enclosed type, -35+50:1°C, capillary prismatic uncoated, special red filling, 260x7-8mm, fully immersed, top with ring

## Additional Information

---

No.	LL-2185
Manufacturer (Brand)	Amarell
VGKL number	161013003
EAN	4058072463229
diameter	8 mm
Temp MAX	50 °C
Temp MIN	-35 °C
Country of origin	Germany
Region of origin	Bavaria
Customs tariff number	90251180
Length	260 mm
Color measuring liquid	red

