



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® refrigerator thermometer

€10.20
plus VAT &
Shipping

Product Images



Description

Simple freezer thermometer for checking the temperature in refrigerators .

With clip for hanging. Should not be missing in any refrigerator.

Colored scale.

Measuring range: -30°C to +30°C

Dimensions: 72 mm Ø

Additional Information

No.	1-2051
Manufacturer (Brand)	neoLab
EAN	4058072033774
diameter	72 mm
Temp MAX	30 °C
Temp MIN	-30 °C
Display (reading, setting)	analog
Counting and measuring devices Features	without mercury
Measuring unit	°C
Division	1 °C
Type thermometer	Analog thermometer

