



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Maxima/Minima thermometer, aluminum case, -30 to +50°C, mercury-free

€17.70
plus VAT &
Shipping

Product Images



Description

Maximum and minimum temperatures always in view!

With large, easy-to-read scale and pushbutton for resetting the thermometer.

Measuring range: -30°C to +50°C

Housing color: gray

Both models are mercury-free.

Dimensions: 220 x 60 mm (L x W)

Additional Information

No.	2-9817
Manufacturer (Brand)	neoLab
EAN	4058072011284
Transport temperature	Room temperature
diameter	3,5 mm
Temp MAX	50 °C
Temp MIN	-30 °C
Wide	60 mm
Height	220 mm
Length	220 mm
Depth	220 mm
Display (reading, setting)	analog
Counting and measuring devices Features	without mercury with maximum/minimum function
Figure	2
Color measuring liquid	red
Measuring unit	°C
Division	1 °C
Type thermometer	Analog thermometer

