



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine thermometer without mercury + protective cover, -35 to + 50°C

€18.10
plus VAT &
Shipping

Product Images



Description

An all-purpose thermometer with special blue liquid. Easy to read. Shatterproof and shapely, stored in metal case, chrome-plated, with metal clip for fastening in any breast pocket or jacket pocket. Two measuring ranges, each divided into 1/1°C

Total length: approx. 12 resp. 15 cm

Additional Information

| | |
|---|--------------------|
| No. | 2-9805 |
| Manufacturer (Brand) | neoLabLine |
| EAN | 4058072011239 |
| Transport temperature | Room temperature |
| diameter | 7 mm |
| Temp MAX | 50 °C |
| Temp MIN | -35 °C |
| Length | 120 mm |
| Display (reading, setting) | analog |
| Counting and measuring devices Features | without mercury |
| Color measuring liquid | blue |
| Measuring unit | °C |
| Division | 1 °C |
| Type thermometer | Analog thermometer |

