



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Amarell pocket thermometer special liquid, -30°C to +50°C, incl. case

**€8.29** As low as

**€5.10**

**plus VAT &  
Shipping**

### Product Images



## Description

---

Multi-purpose thermometer with special red liquid.

Easy to read, wide scale.

-30°C to +50°C in yellow plastic case

-10°C to +100°C in blue plastic case

To minimize the risk of breakage, the perforation on the thermometer head also allows measurements to be taken in the plastic case.

If necessary, the thermometer can simply be removed from the sleeve.

With fastening clip.

Thermometer Ø: approx. 7 mm

Total length: 155 mm

## Additional Information

---

No.	LL-3346
Manufacturer (Brand)	Amarell
old neoLab article no.	2-9809
EAN	4066292052789
Transport temperature	Room temperature
diameter	10 mm
Temp MAX	50 °C
Temp MIN	-30 °C
Length	165 mm
Display (reading, setting)	analog
Counting and measuring devices Features	without mercury
Color measuring liquid	red
Measuring unit	°C
Division	1 °C
Type thermometer	Analog thermometer

