



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Thermometer for indoor and outdoor temperature

€48.10
plus VAT &
Shipping

Product Images



Description

Handy thermometer with large LCD dual display.

- switchable from °C to F
- Memory function for min./max. temperature values
- Measuring range inside: -10°C to +50°C
- Measuring range outside: -50°C to +70°C
- Resolution: 0,1°C
- Accuracy: +/- 1°C

Outdoor sensor with 3 m cable

Power supply: 1 Mignon battery 1.5 Volt

Dimensions: 102 x 110 x 21 mm (L x W x D)

Weight: 133 g (incl. battery)

Additional Information

No.	2-5425
Manufacturer (Brand)	neoLab
EAN	4058072004644
Transport temperature	Room temperature
Color	White
Temp MAX	70 °C
Temp MIN	-50 °C
Wide	21 mm
Height	110 mm
Length	98 mm
Depth	98 mm
Display (reading, setting)	digital
Weight	133 g
for ambient temperature	-10 to 50 °C
Counting and measuring devices Features	with maximum/minimum function with outdoor sensor with alarm function batteries/battery included
Application thermometer	Inside Outside
Resolution °C	0.1 °C
Energy supply	Battery
Functionality	electronic
Measuring unit	°C °F
Measuring accuracy °C	1 °C
Type probe	Air/gas sensor with cable
Type thermometer	Digital thermometer

