



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Amarell Max./Min. thermometer, -50 to +70°C, resolution 0.1°C

**€25.82** As low  
as  
**€20.43**  
plus VAT &  
Shipping

### Product Images



## Description

---

Shapely and functional!

- Temperature memory function
- simultaneous display of indoor and outdoor temperature
- Display illumination
- Switchover °C to °F
- with 3 m outdoor sensor
- with stand and wall mount
- with 1.5 V (AAA) micro battery

Temperature measuring range:

-50°C to +70°C

Resolution: 0.1°C

## Additional Information

No.	LL-0616
Manufacturer (Brand)	Amarell
old neolab article no.	2-5420
EAN	4058072288198
Color	White
Temp MAX	70 °C
Temp MIN	-50 °C
Wide	20 mm
Height	80 mm
Length	55 mm
Depth	55 mm
Display (reading, setting)	digital
Weight	0 g
for ambient temperature	0 to 40 °C
Counting and measuring devices Features	with maximum/minimum function with outdoor sensor batteries/battery included
Resolution °C	0.1 °C
Application thermometer	Inside Outside
Energy supply	Battery
Functionality	electronic
Type probe	Air/gas sensor with cable
Measuring unit	°C °F
Measuring accuracy °C	1 °C
Type thermometer	Digital thermometer

