




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Pocket Penetration Thermometer -40to +200°C

**€29.20**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Pocket thermometer with long probe (125 mm) for simple temperature measurements in liquids, powdery and plastic materials.

- waterproof
- 1-second measuring cycle
- Max-Min-Hold function
- C/F- switchable
- with plastic sleeve + clip
- incl. button cell LR 44

### Technical data

Accuracy:  $\pm 0.8^{\circ}\text{C}$  (-20 to  $100^{\circ}\text{C}$ ), otherwise  $\pm 2^{\circ}\text{C}$

Dimensions: 205 x 20 x 17 mm (L x W x H)

Weight: 56 g

## Additional Information

No.	4-1070
Manufacturer (Brand)	neoLab
EAN	4058072020651
Transport temperature	Room temperature
Color	White
Temp MAX	200 °C
Temp MIN	-40 °C
Wide	17 mm
Height	202 mm
Length	20 mm
Depth	20 mm
Display (reading, setting)	digital
Weight	44 g
Counting and measuring devices Features	waterproof with maximum/minimum function with insertion probe batteries/battery included
Resolution °C	0.1 °C
Application thermometer	Liquids viscous/powdery material
Energy supply	Battery
Functionality	electronic
Type probe	Insertion sensor Submersible sensor
Measuring unit	°C °F
Measuring accuracy °C	0.8 °C
Type thermometer	Digital thermometer

