



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Amarell Digital Pocket Thermometer, -50 to +200°C

€25.46 As low
as
€18.46
plus VAT &
Shipping

Product Images



Description

Measuring range: -50°C to +200°C: 0.1°C

Accuracy: $\pm 2^\circ\text{C}$

Storage of the moment indication possible.

- Stainless steel probe with piercing tip
- complete with battery.
- with on/off switch.
- Battery life up to 5000 hours
- with plastic sleeve, can be used as a handle

Accessories: button cell V 76 PX (SR 44) (order no. 4-1089)

Additional Information

No.	LL-0599
Manufacturer (Brand)	Amarell
old neolab article no.	4-1099
EAN	4058072288181
Color	Black
Temp MAX	200 °C
Temp MIN	-50 °C
Wide	22 mm
Height	147 mm
Length	55 mm
Depth	55 mm
Display (reading, setting)	digital
Weight	0 g
Counting and measuring devices Features	waterproof with insertion sensor batteries/battery included
Application thermometer	Liquids viscous/powdery material
Resolution °C	0.1 °C
Energy supply	Battery
Functionality	electronic
Measuring unit	°C
Measuring accuracy °C	2 °C
Type probe	Insertion sensor Submersible sensor
Type thermometer	Digital thermometer

