



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® pocket thermometer foldable

€32.40
plus VAT &
Shipping

Product Images



Description

Foldable pocket thermometer, waterproof. Ideal for measuring food (according to EN13485), liquids, etc. -40°C to 250°C

Probe tube: 70 mm

Accuracy: $\pm 0.5^{\circ}\text{C}$

Dimensions: 116 x 38 x 20 mm (L x W x H)

Accessories: button cell CR 2032 (order no. 4-1049)

Additional Information

No.	5-2123
Manufacturer (Brand)	neoLab
EAN	4058072024116
Transport temperature	Room temperature
Color	White
Temp MAX	350 °C
Temp MIN	-50 °C
Wide	18 mm
Height	110 mm
Length	35 mm
Depth	35 mm
Display (reading, setting)	digital
Weight	100 g
for ambient temperature	0 to 50 °C
Counting and measuring devices Features	waterproof with insertion sensor batteries/battery included
Resolution °C	0.1 °C
Application thermometer	Liquids viscous/powdery material
Energy supply	Battery
Functionality	electronic
Type probe	Insertion sensor Submersible sensor
Measuring unit	°C °F
Measuring accuracy °C	0.5 °C
Type thermometer	Digital thermometer

