




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Amarell Digital Thermometer with Stainless Steel Probe -35 to +500°C

€70.00 As low
as
€60.56
plus VAT &
Shipping



Description

Extremely fast **digital thermometer**, suitable for industrial use, for obtaining rapid readings even in high temperature ranges. Very fast indicating thermocouple probe made of stainless steel.

Automatic battery test avoids false readings.

Technical features:

- Accuracy: $\pm 1^{\circ}\text{C}$ or 1% of reading
- Battery: Block battery 9V 6F22
- Service life: approx. 2000 h
- Response time: $t_{90\%} = 2$ sec.
- Probe dimensions (l x \varnothing): approx. 100 x 2.5 mm
- Cable: PTFE approx. 400 mm

Additional Information

No.	LL-0621
Manufacturer (Brand)	Amarell
old neoLab article no.	4-1101
EAN	4058072288204
Color	Black
Temp MAX	500 °C
Temp MIN	-35 °C
Wide	60 mm
Height	25 mm
Length	90 mm
Depth	90 mm
Display (reading, setting)	digital
for ambient temperature	0 to 50 °C
Counting and measuring devices Features	with insertion sensor including batteries/battery pack
Resolution °C	1 °C
Application thermometer	Liquids viscous/powdery material
Energy supply	Battery Accumulator
Functionality	electronic
Type probe	Insertion probe
Measuring unit	°C
Measuring accuracy °C	1 °C
Type thermometer	Digital thermometer

