



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Amarell Electronic Digital Thermometer, Maxi-Pen, -50...+200:0.1°C, switchable to

€31.99 As low as
€24.33
plus VAT &
Shipping



Description

Electronic Digital Thermometer, Maxi-Pen, -50...+200:0,1°C, switchable to °F, maximum-minimum-memory, hold-function, insertion probe made of stainless steel 125x3,5mm, water protected

Additional Information

No.	LL-0607
Manufacturer (Brand)	Amarell
VGKL number	681067918
EAN	4058072462437
Temp MAX	200 °C
Temp MIN	-50 °C
Country of origin	China
Customs tariff number	90251900
Warm-up time	4 s
measured size(s)	Temperature

