



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Amarell Electronic Digital Thermometer, Maxi-T, -50...+200:0.1°C, switchable to

**€30.12** As low as  
**€22.97**  
**plus VAT &  
Shipping**

## Product Images

---



## Description

---

Electronic Digital Thermometer, Maxi-T, -50...+200:0,1°C, switchable to °F, maximum–minimum–memory, hold–function, waterproof, stainless steel insertion probe 125x3,5mm

## Additional Information

---

No.	LL-0606
Manufacturer (Brand)	Amarell
VGKL number	681067916
EAN	4058072185299
Temp MAX	200 °C
Temp MIN	-50 °C
Country of origin	China
Customs tariff number	90251900
Warm-up time	4 s
measured size(s)	Temperature

