



## neoLab® Infrared Thermometer

### Product Images

---



### Description

---

For non-contact measurement of surface temperature  
Application in household, hobby and professional areas

Measuring time: 1 second  
Dimensions: 91 x 36 x 14 mm (L x W x H)  
Weight: 50 g  
Operating instructions included

Accessories: Button cell CR 2032 (Order No. 4-1049)

## Additional Information

No.	5-2122
Manufacturer (Brand)	neoLab
EAN	4058072024109
Color	Weiß
Temp MAX	199 °C
Temp MIN	-33 °C
Wide	14 mm
Height	91 mm
Length	36 mm
Depth	36 mm
Display (reading, setting)	digital
Weight	50 g
for ambient temperature	0 bis 50 °C
Counting and measuring devices Features	inklusive Batterien/Akku
Resolution °C	0.5 °C
Application thermometer	Oberflächen
Energy supply	Batterie
Functionality	elektronisch
Type probe	Infrarot
Measuring unit	°C
Measuring accuracy °C	2 °C
Type thermometer	Digital-Thermometer Infrarot-Thermometer

