




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



2mag contact thermometer for magnetic stirrer 7-3080

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Contact thermometer with Pt 1000 probe

Safety and control circuit according to DIN 12878 class 1 and 2, display backlight, LED display for control status, IP65, approx. 150 g, incl. holder for stand rod.

Additional Information

No.	ZU-0058
Manufacturer (Brand)	2mag
old neolab article no.	7-3082
EAN	4058072444464
Transport temperature	Room temperature
Color	light gray
Temp MAX	300 °C
Temp MIN	-50 °C
Wide	22 mm
Height	115 mm
Length	54 mm
Depth	54 mm
Display (reading, setting)	digital
Weight	90 g
Counting and measuring devices Features	waterproof with insertion sensor batteries/battery included
Application thermometer	Liquids
Resolution °C	0.1 °C
Energy supply	Battery
Functionality	electronic
Measuring unit	°C
Measuring accuracy °C	1 °C
Type probe	Immersion sensor with cable
Type thermometer	Digital thermometer

