



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e



Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Amarell laboratory thermometer, DIN 12778, enclosed type, -10/0+100:1°C, capillary

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Laboratory thermometer, DIN 12778, enclosed type, -10/0+100:1°C, capillary prismatic uncoated, red special filling, 305x8-8,5mm, fully immersing, Richter lock with knob, officially calibrated, with calibration certificate with 3 test points

Additional Information

No.	LL-2467
Manufacturer (Brand)	Amarell
VGKL number	101015301
EAN	4058072472429
diameter	8,5 mm
Temp MAX	100 °C
Temp MIN	0 °C
Country of origin	Germany
Region of origin	Bavaria
Customs tariff number	90251180
Length	305 mm
Color measuring liquid	red

