



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657




Amarell pour point thermometer, DIN 12785, enclosed type, -70+50:1°C, capillary

€63.48 As low

as

€48.52

**plus VAT &
Shipping**



Description

Pour point thermometer, DIN 12785, enclosed type, -70+50:1°C, capillary prismatic uncoated, toluene filling, upper part 195x10mm, lower part 165x5mm, immersion depth 180mm, Richter lock without knob, suitable for official certification

Additional Information

No.	LL-2592
Manufacturer (Brand)	Amarell
EAN	4058072474454
Temp MAX	50 °C
Temp MIN	-70 °C
Country of origin	Germany
Region of origin	Bavaria
Customs tariff number	90251180
Color measuring liquid	red

