



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e



Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Amarell thermometer with standard ground joint NS 14.5/23, similar to DIN, enclosed form

€59.11 As low as
€45.60
**plus VAT &
Shipping**

Product Images



Description

Thermometer with standard ground joint NS 14,5/23, similar to DIN, enclosed type, -10+150:1°C, capillary prismatic uncoated, red special filling, upper part 250x11mm, immersion length from upper edge ground joint 75mm, immersion depth from lower edge ground joint 52mm, Richter lock without knob, suitable for official certification

Additional Information

No.	LL-2154
Manufacturer (Brand)	Amarell
VGKL number	161017000
EAN	4058072288280
Temp MAX	150 °C
Temp MIN	-10 °C
Country of origin	Germany
Region of origin	Bavaria
Customs tariff number	90251180
Standard ground core	44695
Color measuring liquid	red

