



Testo 112 Temperature measuring instrument

Product Images

Description

testo 112 Temperature measuring instrument

Additional Information

No.	TO-0907
Manufacturer (Brand)	Testo
EAN	4066292356801
Gross weight	0.32
Manufacturer EAN	4029547003967
Temp MAX	300 °C
Temp MIN	-50 °C
Country of origin	Deutschland
Packing width	0.29000 m
Packing height	0.05000 m
Packing depth	0.11000 m
Wide	64 mm
Customs tariff number	90251900900
Height	182 mm
Depth	40 mm
measured size(s)	Temperatur

