




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Wiggens vacuum controller, measuring accuracy 0.1 mbar

**€2,879.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The precise **Wiggins DVR 480** is a multifunctional **vacuum controller** for connection to pumps or vacuum networks. They offer a wide range of measurement and control in vacuum technology and are ideally suited for applications where an accurate vacuum is essential.

The vacuum is adjustable from 0.1 to 1000 m.

The unit features a bright LED display and convenient on-touch control.

- Suitable for rotary evaporators
- Supports programming of control with up to 5 steps
- All parts in contact with gases are made of PTFE or high-strength ceramics
- High chemical resistance
- RS-232 and analog connection available
- Direct electrical control of the vacuum pump
- Suitable for continuous operation

## Additional Information

---

No.	WI-1545
Manufacturer (Brand)	Wiggins
old neoLab article no.	F-1028
EAN	4058072443221
Transport temperature	Room temperature
Operation via	Membrane keypad
Pressure MAX	1.00 bar abs.
Pressure MIN	0.1000 mbar abs.
Energy supply	Power grid
Measuring unit	mbar
for medium	Air/Gases/Vapors
measured size(s)	Print

