



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® 3-channel multi-timer with large display


Special Price

**€14.50** was

**€22.90**

plus VAT &  
Shipping

### Product Images



## Description

---

3-channel timer with count-down and count-up function

- each channel with separate display
- Clock function, optional 12 or 24 h mode
- Running time each 99 h, 59 min, 59 sec. in second steps
- Start/stop button
- with holding magnet, holding clip and stand

Including 2 x 1.5V lithium battery

## Additional Information

---

No.	2-1999
Manufacturer (Brand)	neoLabLine
EAN	4058072154882
Transport temperature	Room temperature
Color	Grey white orange
Material	Acrylonitrile butadiene styrene (ABS)
Wide	68 mm
Form	rectangular
Height	88.7 mm
Display (reading, setting)	digital
Operation via	Knobs/Buttons
Counting and measuring devices Features	with countup function/stopwatch with countdown function/short timer with addition function batteries/battery included
Number of channels/ timer	3
Number of digits	6
Energy supply	Battery
Functionality	electronic
Type battery	Battery
Type watch	Short timer/ timer. digital

