



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLab® 4 zone timer with note field

**€11.10**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

An intelligent timer with 4 independent timer functions. The white note field offers the possibility to note the different timings of a process. The pen is equipped with a wiping rubber. Batteries type micro AA (2 pcs) included.

Functions:

- maximum time for count down/up:  
24 hours and 59 minutes
- up to 4 different time programs adjustable
- large writing field
- Stand clip and strong magnet on the back
- long beep approx. 5 minutes
- large LCD display (44 x 55 mm)
- Timer housing made of ABS
- Dimensions: 125,5 x 88,5 x 23,5 mm (L x W x D)

Accessories: batteries alkaline micro AAA (order no. 2-3023)

## Additional Information

---

No.	2-2039
Manufacturer (Brand)	neoLabLine
EAN	4058072059385
Transport temperature	Room temperature
Color	yellow
Form	rectangular
Display (reading, setting)	digital
Counting and measuring devices Features	with countup function/stopwatch with countdown function/short timer with addition function batteries/battery included
Number of channels/ timer	4
Functionality	electronic
Spare part	2-2038
Type battery	Battery
Type watch	Short timer/ timer. digital

