



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Timer with direct numerical input

€10.80
plus VAT &
Shipping

Product Images



Description

This digital timer with stopwatch function is especially suitable when exact times must be kept. With the number field you can set any time up to 99 std, 99 min and 99 sec.

- Includes 2 x 1.5 V AAA batteries
- Dimensions: 66 x 20 x 108 mm
- Weight: 50 g

Additional Information

No.	2-4576
Manufacturer (Brand)	neoLab
EAN	4058072159269
Transport temperature	Room temperature
Color	White
Form	rectangular
Display (reading, setting)	digital
Counting and measuring devices Features	with countup function/stopwatch with countdown function/short timer batteries/battery included
Number of channels/ timer	1
Spare part	2-3023
Functionality	electronic
Runtime	99 hrs.
Type battery	1.5V AAA
Type watch	Short timer/ timer. digital

