



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® timer with extra loud alarm

€11.90
plus VAT &
Shipping

Product Images



Description

The neoLab **digital timer with stopwatch function**, fast time input and extra loud alarm is especially needed where there is a relatively high noise level. Whether in production facilities or in laboratories where large devices provide a high noise level, this timer will no longer be overheard. In addition to the alarm, a small light comes on to provide a visual alarm.

- large display
- fast time indication of 10, 5 and 1 min. and 10 s.
- with magnet to attach to refrigerator, workbench,...
- with stand and suspension eye
- Dimensions: 90 x 17 x 90 mm
- Weight: 66 g

Additional Information

No.	2-4574
Manufacturer (Brand)	neoLab
EAN	4058072158392
Transport temperature	Room temperature
Color	White
Form	rectangular
Display (reading, setting)	digital
Counting and measuring devices Features	with countup function/stopwatch with countdown function/short timer batteries/battery included
Number of channels/ timer	1
Functionality	electronic
Spare part	4-1049
Runtime	99 hrs.
Type battery	1.5V AAA
Type watch	Short timer/ timer. digital

