




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Day timer, mechanical

€13.40
plus VAT &
Shipping

Product Images



Description

Mechanical timer with 96 on and off times per day.

- minimum switching interval 15 minutes
- integrated child safety lock
- max. load 16 A, 3600 W
- Rated voltage 230 V, 50 Hz

Additional Information

No.	2-1773
Manufacturer (Brand)	neoLab
EAN	4058072057909
Transport temperature	Room temperature
Display (reading, setting)	analog
Voltage	230 V
Functionality	mechanical electrical
Power	3600 W

