



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Timer up to 10 hours, incl. battery

€13.80
plus VAT &
Shipping

Product Images



Description

Electronic timer and stopwatch with large display.

Direct number entry, up to 10 hours.

With mounting magnet, stand, clip and hanger.

Battery 1.5 V AAA included

Dimensions: 70 x 60 x 25 mm (L x W x D)

Weight: 95 g

Additional Information

No.	2-1399
Manufacturer (Brand)	neoLab
EAN	4058072056506
Transport temperature	Room temperature
Color	Blue White
Form	rectangular
Display (reading, setting)	digital
Counting and measuring devices Features	with countdown function/short timer with addition function batteries/battery included
Number of channels/ timer	1
Functionality	electronic
Type battery	Battery
Type watch	Short timer/ timer. digital

