



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Timer with time program plug for weekly program

Special Price

€37.00 was

€44.70

**plus VAT &
Shipping**

Product Images



Description

Small and fine!

Program plug for daily or weekly programs.

Switch-on and switch-off times adjustable. High temperature resistant and self-extinguishing materials.

Order No. 2-1721

Weekly time program with 2-hour segments.

Simple program setting by toggle segments.

Manual switch for immediate switching on or off of the devices.

Switching status display ON/OFF

Switching capacity: 16 A, 230 V

Rated voltage: 230 V, ± 10%

Additional Information

No.	2-1721
Manufacturer (Brand)	neoLab
EAN	4058072057862
Transport temperature	Room temperature
Display (reading, setting)	analog
Voltage	230 V
Functionality	mechanical electrical
Power	3600 W

