



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hecht table stopwatch 0-60 / 0-100 min., Ø 110 mm

€136.95
plus VAT &
Shipping

Product Images



Description

Large table stopwatch with easy to read dial (Ø 110 mm) and two hands.

- 60 minutes/60 seconds
- or 100/100 minutes
- Time resolution 1 sec. and 1/100 minute
- Hands and number lines luminescent
- Start/ stop/ zero setting
- Addition, summation of time
- Tempo control

Color: white

Dimensions: 170 x 130 x 90 mm (H x W x D)

Weight: 600 g

Accessories: LR6 batteries (order no. 2-3024)

Additional Information

No.	WQ-1628
Manufacturer (Brand)	Hecht
old neoLab article no.	1-1597
EAN	4058072443634
Transport temperature	Room temperature
Color	White
Form	semicircular
Display (reading, setting)	analog
Counting and measuring devices Features	with countup function/stopwatch with addition function
Number of channels/ timer	1
Functionality	electric
Type battery	Battery
Type watch	Short timer. analog

