



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Benchtop Timer, digital

**€82.05**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Table timer with large LCD display (digits 20 mm high).

Normal time display, two timers with separate or common start/stop function, stopwatch, loud alarm tone.

Time display switchable from 999 min., 59 sec. to 9 hrs., 59 min., 59 sec., 1/10 sec.

Incl. mignon battery (AA)

Color: dove blue

Dimensions : 170 x 130 x 80 mm (H x W x D)

## Additional Information

---

|   |   |
|---|---|
| No.                                     | 7-5700  |
| Manufacturer (Brand)                    | neoLab  |
| EAN                                     | 4058072041540   |
| Transport temperature                   | Room temperature  |
| Color                                   | White   |
| Form                                    | semicircular  |
| Display (reading, setting)              | digital   |
| Counting and measuring devices Features | with countup function/stopwatch with countdown function/short timer with addition function batteries/battery included |
| Number of channels/ timer               | 2   |
| Functionality                           | electronic  |
| Type battery                            | Battery   |
| Type watch                              | Short timer/ timer. digital   |

