



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Wall clock with quartz movement, 20 cm Ø, brushed aluminum

€24.30
plus VAT &
Shipping

Product Images



Description

Wall clock with second hand,
case brushed aluminum.

- Diameter 20 cm
- Height 3,5 cm
- large, easy to read numbers
- precision quartz movement accurate to the second
- dust-free glass cover
- Delivery incl. battery 1,5 V AA

Additional Information

No.	2-4575
Manufacturer (Brand)	neoLab
EAN	4058072003753
Transport temperature	Room temperature
diameter	200 mm
Color	silver white
Material	Aluminum
Display (reading, setting)	analog
Counting and measuring devices Features	batteries/battery included
Thickness	35 mm
Functionality	electric
Type watch	Wall clock

