



SI Analytics Lab845Set/BL19pH

Product Images

Description

Laboratory pH meter, measuring parameters pH, mV, ISE, temp., 2-point calibr., micropr., BNC connection. Incl. stand, power supply, BlueLine 19 pH and DIN buffers in ampoules (6 pcs.)

Additional Information

No.	SY-1029
Manufacturer (Brand)	SI Analytics
VGKL number	102706810
EAN	4066292316591
Gross weight	192
Net weight	192.00000 kg
Country of origin	Deutschland
Customs tariff number	90278020

