



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Mettler-Toledo pH meter, FiveEasy F20 meter

**€565.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The new **FiveEasy benchtop instruments** are ideal for laboratory applications and provide high quality pH/mV measurements at the touch of a button.

With their compact design and various data export options, the FiveEasy series instruments are the perfect choice for obtaining reliable measurement data at a reasonable price.

**Without electrode**

## Additional Information

|   |   |
|---|---|
| No.                                     | TD-0168   |
| Manufacturer (Brand)                    | Mettler-Toledo  |
| VGKL number                             | 105006658   |
| EAN                                     | 4066292341395   |
| Transport temperature                   | Room temperature  |
| autoclavable at 121°C                   | No  |
| Color                                   | whitelight gray   |
| Material                                | Acrylonitrile butadiene styrene (ABS)   |
| Wide                                    | 227 mm  |
| Height                                  | 70 mm   |
| Depth                                   | 147 mm  |
| Display (reading, setting)              | Display 4.3" LCD Color  |
| Weight                                  | 0.63 kg   |
| Type interface                          | RS232 USB   |
| for ambient temperature                 | 0 to 40 °C  |
| Counting and measuring devices Features | with internal calibration with integrated data memory with automatic temperature compensation (ATC) with automatic buffer recognition with interface with limit value setting with alarm function |
| Resolution pH                           | 0.01 pH   |
| (Combined) functions                    | pH meter thermometer  |
| Energy supply                           | Power grid  |
| Calibration type                        | at the touch of a button  |
| measured size(s)                        | pH value temperature  |
| pH MAX                                  | 14.0 pH   |
| pH MIN                                  | 0.0 pH  |
| Data transmission                       | USB Cable   |
| Shape connection                        | straight  |
| Manufacturer                            | Mettler-Toledo OnLine   |

|                       |              |
|-----------------------|--------------|
| Measuring accuracy pH | 0.01 pH      |
| Voltage MAX           | -1.999900 mV |
| Voltage MIN           | 1.999900 mV  |
| Type plug             | BNC          |

