



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Mettler-Toledo conductivity meter FiveGo, F3 meter

€423.00
plus VAT &
Shipping

Product Images



Description

The **FiveGo™ portable instruments** are ideal for mobile applications in the laboratory and provide high quality measurements for conductivity with just one push of a button.

The rugged and waterproof design makes the FiveGo™ instruments perfect for obtaining reliable measurement data at a reasonable price, even under demanding conditions.

Without electrode

Additional Information

No.	TD-0171
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105006883
EAN	4066292341425
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	whitelight gray
Material	Acrylonitrile butadiene styrene (ABS)
Temp MAX	100 °C
Temp MIN	0 °C
Wide	77 mm
Height	188 mm
Depth	33 mm
Weight	0.26 kg
for ambient temperature	0 to 40 °C
Counting and measuring devices Features	waterproof with internal calibration with integrated data memory with limit setting with alarm function batteries/battery included
Resolution °C	0.1 °C
(Combined) functions	Conductivity meter Thermometer
Energy supply	Battery Battery
Calibration type	at the touch of a button
measured size(s)	Conductivity TDS in g/l Temperature
Data transmission	Cable
Shape connection	straight
Manufacturer	Mettler-Toledo OnLine
Measuring accuracy °C	0.5 °C
Voltage MAX	1999.000000 mV
Voltage MIN	-1999.000000 mV

