



Mettler-Toledo SevenCompact S210 kit

Product Images

Description

High-performance routine instrument kit for pH, mV, rel. mV, and temperature measurement: - Large color display (4.3 inch) with u-Focus™: easily switch display all/only the important information - Routine or expert mode - Full text menu navigation in 10 languages - ID management: User/sample/sensor ID/full GxP support - Automatic 1-/2-/3-/4-/5-point calibration with 8 predefined and user-defined buffer groups - Calibration timer / reminder - Automatic, manual or time-controlled end point acquisition - adjustable end point criteria - real-time clock and date - internal measured value memory for 1000 data sets - 2 USB interfaces for data export via USB stick /PC, connection of USB scanner(s) or keyboard possible - control of a stirrer (optional accessory) - U-Place™: innovative, absolutely vertical running electrode holder => avoids knocking over of vessels and electrode breakage - IP54 dust- and. Moisture-proof - and much more. Scope of delivery: Instrument, electrode InLab Expert Pro ISM holder, set of buffer solutions, operating instructions, declaration of conformity and test certificate Technical data: - pH measuring range: -2.000 to 20.000, resolution 0.001 /0.01 / 0.1

Additional Information

No.	TD-0109
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105009026
EAN	4066292340824
Gross weight	3
Country of origin	China
Customs tariff number	90278020

