



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® XS® pH 1 Tester ECOPACK

€71.00
plus VAT &
Shipping

Product Images



Description

Portable pH meter waterproof IP67

PH 1 meter series with fixed and non-replaceable electrode

PH 5 meter series with replaceable electrode

- Resolution 0.1 pH with microprocessor and fixed sensor
- Calibration with key for 1 or 2 points (selection of pH 4, 7, 10)
- Large LCD display with indication of buffers used for calibration and measurement stability
- Automatic temperature compensation
- Automatic shutdown
- Messages for self-diagnosis: batteries, errors in calibration
- Power supply 2 x 1.5V AAA batteries, also compatible with rechargeable batteries
- The ECO PACK includes: additional glass for calibration, carrying strap, quick reference guide and 2 x 1.5V AAA batteries.

Measuring range: 0.0 ... 14.0 pH (pH 0.1). Working temperature: 0..50 ° C

Additional Information

No.	4-3070
Manufacturer (Brand)	neoLab
EAN	4058072190552
Transport temperature	Room temperature
diameter	0,1 mm
Color	green
Length	200 mm
Display (reading, setting)	graphic display
Weight	0.11 kg
Resolution pH	0.01 pH
(Combined) functions	pH meter
Energy supply	Battery
pH MAX	14.0 pH
pH MIN	0.0 pH
Measuring unit	pH
Measuring accuracy pH	0.10 pH
Type battery	1.5V AAA
Type sensor	XS Standard for tester

