



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® pH Check, pH 0 to 14, 0 ° to 50 °C,
batteries included**

**€53.30
plus VAT &
Shipping**

Product Images



Description

Compact pH tester for a range from pH 0 to 14.00 and a temperature range from 0 °C to 50 °C

- Ideally suited for measurement in narrow vessels such as test tubes
- IP65 splash-proof
- With simple one-touch calibration function and hold function
- Split display with automatic temperature compensation

Additional Information

No.	4-1062
Manufacturer (Brand)	neoLab
EAN	4058072157357
Transport temperature	Room temperature
Color	Blue White
Temp MAX	50 °C
Temp MIN	0 °C
Wide	20 mm
Height	225 mm
Length	35 mm
Depth	35 mm
Display (reading, setting)	digital
Weight	69 g
Counting and measuring devices Features	waterproof with internal calibration with certification batteries/battery included
Resolution pH	0.01 pH
Resolution °C	0.1 °C
(Combined) functions	pH meter
Energy supply	Battery
Calibration type	at the touch of a button
measured size(s)	pH value
pH MAX	14.0 pH
pH MIN	0.0 pH
Measuring unit	°C °F pH
Measuring accuracy pH	0.10 pH
Measuring accuracy °C	1 °C
Type certificate	Declaration of conformity (standardization)

