




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Thermo Scientific pH pocket tester, 1.0 to 15.0, with ATC

**€199.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

These pH testers are designed for a long service life with high measuring accuracy and low operating costs. The double contact electrodes are ideal for use in dirty water or in solutions with heavy metals or other contaminants, also because the increased volume of the reference gel extends the service life of the electrode.

**OrderNo. 8-1302**

pHTestr20 with ATC

**Measuring range: 1.0 to 15.0 pH**

**Accuracy:  $\pm 0.01$  pH**

Accessories: Spare electrode (Order No. 8-1304)

## Additional Information

No.	LD-4918
Manufacturer (Brand)	Thermo Scientific
old neolab article no.	8-1302
EAN	4066292041950
Color	Yellow White
Wide	40 mm
Height	40 mm
Length	165 mm
Depth	165 mm
Display (reading, setting)	digital
Weight	125 g
for ambient temperature	0 to 50 °C
Counting and measuring devices Features	waterproof with automatic temperature compensation (ATC) batteries/battery included
(Combined) functions	pH meter
Energy supply	Battery
Calibration type	automatic
measured size(s)	pH value
pH MAX	15.0 pH
pH MIN	1.0 pH
Spare part	LD-4911
Measuring unit	pH
Measuring accuracy pH	0.01 pH

