



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific pH pocket tester, 1.0 to 15.0, with ATC & temperature display

€214.00
plus VAT &
Shipping

Product Images

Description

These pH testers are designed for a long service life with high measurement accuracy and low operating costs. The double contact electrodes are ideal for use in dirty water or in solutions with heavy metals or other contaminants, also because the increased volume of the reference gel extends the service life of the electrode.

Order No. 8-1303

pHTestr 30 with ATC and temperature display

Measuring range: 1.0 to 15.0 pH

Accuracy: ± 0.01 pH

Accessories: Spare electrode (Order No. 8-1304)

Additional Information

No.	LD-4919
Manufacturer (Brand)	Thermo Scientific
old neolab article no.	8-1303
EAN	4066292041967
Color	Yellow White
Temp MAX	50 °C
Temp MIN	0 °C
Wide	40 mm
Height	40 mm
Length	165 mm
Depth	165 mm
Display (reading, setting)	digital
Weight	125 g
for ambient temperature	0 to 50 °C
Counting and measuring devices Features	waterproof with automatic temperature compensation (ATC) batteries/battery included
Resolution pH	0.01 pH
Resolution °C	0.1 °C
(Combined) functions	pH meter thermometer
Energy supply	Battery
Calibration type	automatic
measured size(s)	pH value temperature
pH MAX	15.0 pH
pH MIN	1.0 pH
Spare part	LD-4911
Measuring unit	pH °C °F
Measuring accuracy pH	0.01 pH

