



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton Micro Tissue Hand Homogenizer Tenbroek 1 ml, 2 pcs/pack

€231.28
plus VAT &
Shipping

Product Images

Description

A complete selection of Wheaton tissue homogenizers as a set of 7 different models.
In plastic case with shovel insert for safe storage.

The set contains 1 each

- 0.1 ml glass pestle (9-0905)
- 0.2 ml glass pistil with Teflon® ring (9-0906)
- 0.5 ml Teflon® pistil, with screw cap (9-0900)
- 1 ml Tenbroek, glass pestle (9-0971)
- 1 ml Dounce, glass pestle (9-0972)
- 1 ml Teflon® pestle, tapered (9-0950)
- 2 ml Potter-Elvehjem Teflon® pestle (9-0973)

The glass-to-glass surface homogenizes connective cell tissues.

The Teflon® units are particularly suitable for soft tissues.

No. 9-0900 is suitable for extraction, separation or purification with a 13 mm screw cap.

Additional Information

No.	WD-1284
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0971
EAN	4058072427351
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume	1 ml
Functionality	mechanical
Number of pistils	1
Material rod	Borosilicate glass
DM Mortar	11 mm
DM Rod	10 mm
DM Stamp	10.91 mm
Shape stamp	around
Length mortar	48 mm
Length pestle	140 mm
Material mortar	Borosilicate glass
Material stamp	Borosilicate glass
Model Tissue Homogenizer Set	Tenbroeck (2x)
Type homogenizer	Tenbroeck
Wall distance pestle 1	0.09 to 0.16 mm

