



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton Micro Tissue Hand Homogenizer Potter Elvehjem 2 ml, 2pcs/pack

€315.80
plus VAT &
Shipping

Product Images

Description

A complete selection of Wheaton tissue homogenizers as a set of 7 different models. In plastic case with shovel insert for safe storage.

The set contains 1 each

- 0.1 ml glass pestle (9-0905)
- 0.2 ml glass pistil with Teflon® ring (9-0906)
- 0.5 ml Teflon® pistil, with screw cap (9-0900)
- 1 ml Tenbroek, glass pestle (9-0971)
- 1 ml Dounce, glass pestle (9-0972)
- 1 ml Teflon® pestle, tapered (9-0950)
- 2 ml Potter-Elvehjem Teflon® pestle (9-0973)

The glass-to-glass surface homogenizes connective cell tissues.

The Teflon® units are particularly suitable for soft tissues.

No. 9-0900 is suitable for extraction, separation or purification with a 13 mm screw cap.

Additional Information

No.	WD-1329
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0973
EAN	4058072427559
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume	2 ml
Number of pistils	1
DM Mortar	11 mm
DM Rod	6.3 mm
DM Stamp	10.9 mm
Shape stamp	around
Functionality	mechanical
Length mortar	45 mm
Length pestle	203 mm
Material mortar	Borosilicate glass
Material rod	Stainless steel
Material stamp	Polytetrafluoroethylene (PTFE. Teflon)
Model Tissue Homogenizer Set	Potter Elvehjem (2x)
Type homogenizer	Potter-Elvehjem
Wall distance pestle 1	0.1 to 0.15 mm

