



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton micro tissue hand homogenizer cpl. set in case, 7 pcs./set

€818.53
plus VAT &
Shipping

Product Images



Description

A complete selection of Wheaton tissue homogenizers as a set of 7 different models. In plastic case with shovel insert for safe storage.

The set contains 1 each

- 0.1 ml glass pestle (9-0905)
- 0.2 ml glass pistil with Teflon® ring (9-0906)
- 0.5 ml Teflon® pistil, with screw cap (9-0900)
- 1 ml Tenbroek, glass pestle (9-0971)
- 1 ml Dounce, glass pestle (9-0972)
- 1 ml Teflon® pestle, tapered (9-0950)
- 2 ml Potter-Elvehjem Teflon® pestle (9-0973)

The glass-to-glass surface homogenizes connective cell tissues.

The Teflon® units are particularly suitable for soft tissues.

No. 9-0900 is suitable for extraction, separation or purification with a 13 mm screw cap.

Additional Information

No.	WD-1369
Manufacturer (Brand)	Wheaton
old neolab article no.	9-0970
EAN	4058072427757
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Shape stamp	round conical. pointed
Functionality	mechanical
Material mortar	Borosilicate glass
Material rod	Borosilicate glass stainless steel
Material stamp	Borosilicate glass Polytetrafluoroethylene (PTFE. Teflon)
Model Tissue Homogenizer Set	Complete set in case
Type homogenizer	Dounce Potter-Elvehjem Tenbroeck Handheld homogenizer. manual

