



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton Micro Tissue Hand Homogenizer, tapered, glass/glass, 5 ml

€177.78
plus VAT &
Shipping

Product Images



Description

Pointed micro hand homogenizer for gradual size reduction.
Coarse size is initially reduced in the upper part, fine in the lower part.

Additional Information

| | |
|------------------------|----------------------------|
| No. | WD-1353 |
| Manufacturer (Brand) | Wheaton |
| old neoLab article no. | 9-0937 |
| EAN | 4058072427634 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Volume | 5 ml |
| Number of pistils | 1 |
| DM Mortar | 13 mm |
| DM Rod | 6.3 mm |
| DM Stamp | 12.9 mm |
| Shape stamp | conical. pointed |
| Functionality | mechanical |
| Length mortar | 93 mm |
| Length pestle | 220 mm |
| Material mortar | Borosilicate glass |
| Material rod | Borosilicate glass |
| Material stamp | Borosilicate glass |
| Type homogenizer | Manual homogenizer. manual |
| Wall distance pestle 1 | 0.1 to 0.15 mm |

