



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Wheaton Micro Tissue Hand Homogenizer Tapered, 3 ml

**€157.75**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Tapered micro hand homogenizer, complete

For gradual comminution especially for soft tissues.

Coarse comminution takes place first in the upper part, fine in the lower part.

## Additional Information

No.	WD-1362
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0953
EAN	4058072427696
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume	3 ml
Hand homogenizers and pestles features	can be clamped into stirrer pestle with ring/distance piece
Functionality	mechanical
Number of pistils	1
Material rod	Stainless steel
DM Mortar	11 mm
DM Rod	6.3 mm
DM Stamp	10.9 mm
Spare part	WD-1361 WD-1360
Shape stamp	conical. pointed
Length mortar	86 mm
Length pestle	220 mm
Material mortar	Borosilicate glass
Material stamp	Polytetrafluoroethylene (PTFE. Teflon)
Type homogenizer	Manual homogenizer. manual
Wall distance pestle 1	0.1 to 0.15 mm

